	Notice d'utilisation - Instruction for use Produktinformation - Instrucciones de uso	Réf : IFU- 0340
	MILIEUX PRETS A L'EMPLOI/ READY-TO-USES MEDIA MUELLER HINTON RAPIDE MHR-SIR® - REF. 020502	VERSION 11 Validé : 30/10/25

FR : NOM DU/DES PRODUIT(S) EN: PRODUCT(S) NAME ES: DENOMINACIÓN DEL PRODUCTO(S) IT : NOME DEL PRODOTTO DE: PRODUKT NAME	MUELLER HINTON RAPIDE MHR-SIR®
FR : REFERENCE(S) COMMERCIALE(S) EN: SALES REFERENCE(S) ES: REFERENCIAS COMERCIALES IT: RIFERIMENTO COMMERCIALE DE : GESCHÄFTSREFERENZ	020502

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ENGLISH

INTENDED USE

Rapid culture media designed for study of micro-organism antimicrobial susceptibility testing.

INTRODUCTION

Mueller Hinton medium ⁽¹⁾, standardized according to W.H.O. guidelines ⁽²⁾, is recommended for antimicrobial susceptibility testing by agar diffusion (Kirby-Bauer method ^(2, 3)) in order to determine susceptibilities of fast growth aerobic microorganisms to antimicrobial agents.

PRODUCT DESCRIPTION

MHR-SIR® medium is composed of a standardized Mueller Hinton agar ^(1,2), supplemented by specific complexes in order to obtain results after 5 to 7 hours of incubation, for fast growing bacteria, such as Enterobacteriaceae and Staphylococcus, under conditions of use and reading specified below.

MHR-SIR® medium can be used only with SIRscan® 2000 automatic and SIRscan® ORION.

COMPOSITION

Approximate formula*

Beef extract	2.0 g
Caseine acid hydrolysate.....	17.5 g
Corn starch.....	1.5 g
Growth factor complex.....	5.5 g
Contrast complex.....	0.05 g
Agar	17.0 g
Purified water.....	1 l

Final pH: 7.2 - 7.4

** This theoretical formula can be adjusted and / or supplemented according to the imposed performance criteria.*

PACKAGING

Ready-to-use media Rapid Mueller Hinton (MHR-SIR®), available in box of 20 square Petri dishes (120x120 mm)

Reference 020502

STORAGE

MHR-SIR® media must be stored between 4 °C and 15 °C, in their box, from receipt and before the expiry date.

Once the bag has been opened, the media must be stored upside down in accordance with CASFM /EUCAST storage conditions.

Expiry date and lot number are indicated on the box label.

	Notice d'utilisation - Instruction for use Produktinformation - Instrucciones de uso	Réf : IFU- 0340
	MILIEUX PRETS A L'EMPLOI/ READY-TO-USE MEDIA MUELLER HINTON RAPIDE MHR-SIR® - REF. 020502	VERSION 11 Validé : 30/10/25

METHOD

Diffusion method for determining the susceptibility of micro-organisms to antibiotics (Kirby-Bauer ^(2,3) method)

PROCEDURE

Equipment provided:

- MHR-SIR® media (ref. 020502)

Equipment required but not supplied:

- SIRscan® prep kit (ref. 041003 – 041004)
- Discs impregnated with antibiotics, SIRscan® Discs range for example
- Antibiotic disc dispensers for round media Ø 90mm (ref. 800226) and for 120x120mm square media (ref. 800256) or forceps
- PreLUD® instrument
- SIRscan® 2000 automatic analyzer or SIRscan® ORION
- Other culture media, reagents and other laboratory equipment required

Inoculum preparation

From colonies isolated on culture medium, prepare a suspension using the Inoclic® system, in tube of physiological water supplied in the SIRscan® prep kit.

Protocol for using Inoclic® systems is specified in its product sheet ⁽⁴⁾.

Inoculation (by swabbing method)

Inoculation can be done manually or via a PreLUD® analyzer equipped with the swab module (or antibiogram module).

- Manual inoculation

Dip a SIRscan® swab in the bacterial suspension and remove excess fluid by rotating the swab on the tube walls.

Streak the entire surface of the MHR-SIR® medium tri-directionally.

Inoculation from a standardized inoculum is detailed in CA-SFM ⁽⁵⁾, EUCAST ⁽⁶⁾ and CLSI ⁽⁷⁾ guidelines.

- Inoculation via PreLUD®

PreLUD® system is defined with a validated method for the preparation of antibiograms in agar diffusion.

Refer to the PreLUD® manual.

Dispensing of antibiotic discs

Discs dispensing can be done manually or automatically using a PreLUD® equipped with the antibiogram module (swab module and discs dispensing).

- Manual dispensing

Apply antibiotic discs with a suitable dispenser or manually with a clamp.

Procedures for antibiotic discs dispensing are defined in the recommendations of CA-SFM ⁽⁵⁾, EUCAST ⁽⁶⁾ and CLSI ⁽⁷⁾.

- Dispensing via PreLUD®

	Notice d'utilisation - Instruction for use Produktinformation - Instrucciones de uso	Réf : IFU- 0340 VERSION 11
	MILIEUX PRETS A L'EMPLOI/ READY-TO-USES MEDIA MUELLER HINTON RAPIDE MHR-SIR® - REF. 020502	Validé : 30/10/25

PreLUD® system with the disc dispensing module is configured and validated for dispensing SIRscan® discs on round or square plates, inoculated by it.

Refer to the PreLUD® manual.

Incubation

Incubate antibiograms in SIRscan® 2000 automatic or SIRscan® ORION incubator by clearly identifying their nature, MHR-SIR® media.

Incubation of antibiograms on MHR-SIR® media lasts 5 to 7 hours. After this period of incubation, SIRscan® 2000 automatic or SIRscan® ORION reads regularly the plates until the final result is achieved.

In case of micro-organism weak growth, particular results or a concern about laboratory organization making fast reading impossible, incubation can be prolonged until 16-24 hours without significant change of antibiotic susceptibility testing results.

Incubation specificities, according to bacterial species considered, are described in the recommendations of CA-SFM ⁽⁵⁾, EUCAST ⁽⁶⁾ and CLSI ⁽⁷⁾.

Inhibition zone measurement

After appropriated incubation, measurement of inhibition zone diameters around antibiotic discs (with millimeter accuracy) is performed using the SIRscan® 2000 automatic or SIRscan® ORION configured for reading MHR-SIR® media. This setting must be performed by an authorized technician.

Recommended protocol and specific recommendations for reading inhibition zone diameters are detailed in the recommendations of CA-SFM ⁽⁵⁾, EUCAST ⁽⁶⁾ and CLSI ⁽⁷⁾.

RESULTS INTERPRETATION

Interpretation of inhibition zone diameters, for antimicrobial susceptibility tests, is carried out by referring to CA-SFM ⁽⁵⁾, EUCAST ⁽⁸⁾ and CLSI ⁽⁹⁾ recommendations.

Interpretation of obtained inhibition zone diameters, conducted by the SIRWEB® / SIRxpert® softwares, determines clinical categorization of a given microorganism to a specific antibiotic.

PERFORMANCES AND QUALITY CONTROL

Performances

Antibiotic susceptibility results obtained on MHR-SIR® media are 96.8% concordant with results obtained for antibiograms of Enterobacteriaceae performed on conventional Mueller Hinton media ⁽¹⁰⁾.

Quality Control

- Appearance of ready-to-use media MHR-SIR® is a light amber color agar.
- MHR-SIR® media respond to the following cultural performances:

Strains	Expected results
<i>Escherichia coli</i> ATCC® 25922™	Growth
<i>Staphylococcus aureus</i> ATCC® 25923™	Growth

	Notice d'utilisation - Instruction for use Produktinformation - Instrucciones de uso	Réf : IFU- 0340 VERSION 11
	MILIEUX PRETS A L'EMPLOI/ READY-TO-USE MEDIA MUELLER HINTON RAPIDE MHR-SIR® - REF. 020502	Validé : 30/10/25

- Quality control of MHR-SIR® medium performances is carried out using particular couples of strains / antibiotic discs for which inhibition diameters must comply with breakpoints defined in databases delivered with SIRWEB® / SIRxpert® softwares.

PRECAUTIONS OF USE

- In vitro diagnostic reagents for professional use only.
- Single use reagents.
- Do not use reagents after the expiration date.
- Do not use culture media showing signs of microbial contamination or any other sign of deterioration.
- Before use, allow media and tubes to be reheated for inoculum preparation at ambient temperature.
- If excess moisture is observed on the surface, it is advisable to dry the cans for 10 to 30 minutes in an oven at 35°C before use.
- SIRscan® prep kit can be used for bacterial inoculum preparation and antibiogram inoculation.
- Antibiograms must be read by a SIRscan® 2000 automatic or SIRscan® ORION analyzer, selecting the MHR-SIR® medium as media type. The system must have been previously set by an authorized technician.
- Follow at all times current techniques and precautions regarding protection against biological hazards.
- After use, remove used and contaminated reagents according to the internal procedures and in accordance with the legislation in force.

LIMITATIONS OF THE PROCEDURE

- It is recommended to perform antibiograms from fresh culture (less than 24 hours).
- An important number of factors, such as handling time at each protocol step, could affect final results. It is imperative to follow protocols described in SIRscan® reagents instructions for use (including this document), in SIRscan® 2000 automatic or SIRscan® ORION user manual and in CA-SFM ⁽⁵⁾, EUCAST ^(6,8) and CLSI ^(7,9) recommendations.

	Notice d'utilisation - Instruction for use Produktinformation - Instrucciones de uso	Réf : IFU- 0340 VERSION 11
	MILIEUX PRETS A L'EMPLOI/ READY-TO-USE MEDIA MUELLER HINTON RAPIDE MHR-SIR® - REF. 020502	Validé : 30/10/25

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REFERENZEN

- 1 - Mueller J.H. & Hinton J. (1941) A Protein-Free Medium for Primary Isolation of the Gonococcus and Meningococcus. Proc. Soc. Exp. Biol. Med. 48(1): 330-333.
- 2 - World Health Organization (1982) W.H.O. Expert Committee on Biological Standardization. Thirty-second Report. Technical Report Series 673.

	Notice d'utilisation - Instruction for use Produktinformation - Instrucciones de uso	Réf : IFU- 0340 VERSION 11
	MILIEUX PRETS A L'EMPLOI/ READY-TO-USES MEDIA MUELLER HINTON RAPIDE MHR-SIR® - REF. 020502	Validé : 30/10/25

3 - Bauer A.W., Kirby W.M.M., Sherris J.C. & Turck M. (1966) Antibiotic susceptibility testing by a standardized single disk method. Am. J. Clin. Pathol. 45(4):493-496.

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8 - European Committee on Antimicrobial Susceptibility Testing (2018) Breakpoint tables for interpretation of MICs and zone diameters.

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